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
Resonance Ltd.	Operating Manual	Page 9 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

Appendix 1. SPECTROMETER SPECIFICATIONS

Date: Jan. 19, 2004

Initial:

Identification:

Spectrometer Model/Serial Number: RMDI /040121

Entrance Slit:

Type	Fixed
Height/Width	1 mm/0.025mm

Exit Detector:

Type	Linear CCD
Pixel	200 x 14 microns
Array Width	28 mm
Spectral coverage	295 to 437 nm
Spectral coverage with UV filter	295 to 385 nm

Grating:

Grooves/mm	2,400/mm
UV type	
Resolution/System	< 0.4 nm first order

Appendix 2: SPECTROMETER CALIBRATION


Date: Jan. 19,2004

Initial:

Equipment:

Spectrometer Model/Serial Number RMDI /

..... Fixed 25 micron slit

Resonance Ltd.	Operating Manual	Page 11 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

Appendix 3: Script File Listing

```
'SCRIPT FILE NAME RMDI Script 6.bas
' The following script was based on the excelinterface.bas which was script file in Ocean
Optics Platinum Software package
' Script files are written in The Sax Basic Language which provides the core language
definition.
' Sax Basic is Visual Basic for Applications(TM) compatible.
' Sax Basic and Sax Basic Editor are products of Polar Engineering And Consulting

' BEORE RUNNING THIS SCRIPT, YOU MUST ADD A REFERENCE TO THE
' MICROSOFT EXCEL TYPE LIBRARY. YOU MUST SELECT EDIT | REFERENCES.
and check off the reference
' compatible with the version of Excel run on the host computer e g. MICROSOFT EXCEL
XX OBJECT LIBRARY (YY)
'
' The resonanceexcelinterface1.bas script used the MICROSOFT EXCEL 9.0 OBJECT
LIBRARY(1.3)
'
' ALSO, MAKE SURE EXCEL IS OPEN PRIOR TO RUNNING THE MACRO AND
THE EXCELTEST.XLS WORKBOOK (or Worksheet in earlier versions of Excel)
' IS OPEN IN EXCEL BEFORE YOU EXECUTE THIS MACRO

' LOG OF Changes AND ADDITIONS TO THE EXCEL INTERFACE SOFTWARE (by
Bill Morrow bmorrow@resonance.on.ca on June 2,03)
' 1. change program name to resonanceexcelinterface
' 2. set references as above xx = 9.0 and (1.3) part of reference. This process is self
explanatory BUT
' it is vitally important not to use two excel references at once. (e g. Excel 5.0 and
9.0 references
' simultaneously
' 3. Add in integrate module which integrates over wavelength ranges.
' 3a. Use Simpsons rule integration (A MODULE IN OO). Normalize each signal to
1 nanometer by dividing by wl2-wl1
' 4. Set up five wavelengths which includes two background wavelengths at 281 and
292 nm (to estimate of scattered light)
' and three wavelengths around an SO2 absorption band at 309.6 nm
' 5. Calculate band strength by averaging signal on either side and taking ratio with
signal on line
' 6. Use a linear extrapolation of the signal at 292 (based on the slope between 281
and 292) to estimate scattered light
' 7. Remove spectra recording in excel
```

Resonance Ltd.	Operating Manual	Page 12 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

```
'=====
'=====
' DO NOT MODIFY THE FOLLOWING BLOCK OF CODE
Option Explicit
'#uses "classes\menuverify.bas"
Global ooi As OOIBase32Platinum
Global Ready As Integer
Global Const INITDONE=66
Global InEvent As Boolean
' END NONMODIFYABLE BLOCK
'=====
'=====
```

```
'#####
'#####
' Insert global variables here
'#####
'#####
```

```
' declare all excel objects

Global oXL As Excel.Application
Global oWB As Excel.Workbook
Global oSheet As Excel.Worksheet

' define global arrays for wavelength and intensity

Global i As Integer
Global j As Integer
Global jmax As Integer
Global NWL1 As Integer
Global NWL2 As Integer
```

```
'Global PixCount As Integer
```

```
Global wl(2048) As Double
Global data(2048) As Double
```

```
Global Spect_StWL As Double
Global Spect_EnWL As Double
```

```
Global Spect_Stpix As Integer
Global Spect_Enpix As Integer
```

Resonance Ltd.	Operating Manual	Page 13 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

Global Spect_Array(2048) As Integer

Global chan As SPECCHANNEL

Global WLO As Double 'Wavelength offset parameter
 Global WLM As Double 'Wavelength multiplier parameter
 Global CWL1(15) As Double 'Center Wavelength array for SO2 slit set (nanometre units)
 Global CWL2(15) As Double 'Center Wavelength array for Ozone slit set (nanometre units)
 Global SHW1(15) As Double 'Slit HALF width array for SO2 slit set (nanometre units)
 Global SHW2(15) As Double 'Slit HALF width array for Ozone slit set (nanometre units)

Global Stpix1(15) As Integer
 Global Enpix1(15) As Integer
 Global Stpix2(15) As Integer
 Global Enpix2(15) As Integer

Global BaseRat309 As Double
 Global BaseRat311 As Double
 Global BaseRat313 As Double
 Global StartTimer As Double

Public Sub EventGenerator_MacroBegin(ByVal modname As String)

```
'=====
'=====
' DO NOT MODIFY THE FOLLOWING BLOCK OF CODE
'   If Ready=INITDONE Then Exit Sub
'   InEvent=False
'   Set ooi=New OOIBase32Platinum
'   END NONMODIFYABLE BLOCK
'=====
'=====
```

```
'#####
'#####
' Insert your start-up code here
'#####
'#####
```

```
' get the wavelength array
ooi.SpectralData.GetWavelengths(chan,wl(0))
```

Resonance Ltd.	Operating Manual	Page 14 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

j = 0

jmax = 10000

' set starting and ending wavelengths

'Compressed Spectrum Mask

Spect_StWL = 308
Spect_EnWL= 356.03

'SO2 MASK

WLO = 0	'Wave length offset can be
used to correct for temp drift	
WLM = 1	'Wave length span can be used
to correct for temp drift	
NWL1 =9	'Number of slits
CWL1(1) = 270*WLM+WLO	'Center wavelength 1 background reference
SHW1(1) = .25*WLM	'Slit 1 half width corrected for wave
length span	
CWL1(2) = 280*WLM+WLO	'Center wavelength 2 background reference
SHW1(2) = .25*WLM	'Slit 2 half width
CWL1(3) = 309.77*WLM+WLO	'Center wavelength 3 low SO2 absorption
SHW1(3) = 0.25*WLM	'Slit 3 half ...
CWL1(4) = 310.75*WLM+WLO	'Center wavelength 4 high SO2 absorption
SHW1(4) = 0.25*WLM	'Slit 4 half width ...
CWL1(5) = 312.1*WLM+WLO	'Center wavelength 5 low SO2 absorption
SHW1(5) = 0.25*WLM	'Slit 5 half width ...
CWL1(6) = 313.08*WLM+WLO	'Center wavelength 6 high SO2 absorption
SHW1(6) = 0.25*WLM	'Slit 6 half width ...
CWL1(7) = 314.16*WLM+WLO	'Center wavelength 7 low SO2
absorption	
SHW1(7) = 0.25*WLM	'Slit 7 half width ...

Resonance Ltd.	Operating Manual	Page 15 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

absorption

CWL1(8) = 315.43*WLM+WLO 'Center wavelength 8 high SO2

SHW1(8) = 0.25*WLM 'Slit 8 half width ...

CWL1(9) = 316.52*WLM+WLO 'Center wavelength 9 low SO2 absorption

SHW1(9) = 0.25*WLM 'Slit 9 half width ...

BaseRat309 = 1.22

BaseRat311 = 1.01

BaseRat313 = 1.01

'OZONE MASK

NWL2 = 5

CWL2(1) = 306.3*WLM+WLO

SHW2(1) = 0.25*WLM 'Slit 1 half width corr for wave length span

CWL2(2) = 310.0*WLM+WLO

SHW2(2) = 0.25*WLM 'Slit 2 half width...

CWL2(3) = 313.5*WLM+WLO

SHW2(3) = 0.25*WLM 'Slit 3 half width...

CWL2(4) = 316.8*WLM+WLO

SHW2(4) = 0.25*WLM 'Slit 4 half width...

CWL2(5) = 320*WLM+WLO

SHW2(5) = 0.25*WLM 'Slit 5 half width...

' set spectrometer channel
chan=chMaster

' get starting and ending pixel from wavelengths

'For compressed spectrum

Spect_Stpix = ooi.Spectrometer.wavelength.GetPixel(chan,Spect_StWL)
Spect_Enpix = ooi.Spectrometer.wavelength.GetPixel(chan,Spect_EnWL)

'For SO2 mask

For i = 1 To NWL1

Resonance Ltd.	Operating Manual	Page 16 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

Stpix1(i) = ooi.Spectrometer.wavelength.GetPixel(chan,CWL1(i)-SHW1(1))
 Enpix1(i) = ooi.Spectrometer.wavelength.GetPixel(chan,CWL1(i)+SHW1(1))

Next i

'For O3 mask

For i = 1 To NWL2

Stpix2(i) = ooi.Spectrometer.wavelength.GetPixel(chan,CWL2(i)-SHW2(i))
 Enpix2(i) = ooi.Spectrometer.wavelength.GetPixel(chan,CWL2(i)+SHW2(i))

Next i

'Configure OO Sheet

' create new grid window

' If ooi.SetGridByName("Integration") = False Then
 ' ooi.Windows.NewGridWindow(4,30,1,"Integration")
 ' Call ooi.SetGridActive()
 ' End If

' name the cells

' ooi.Grid.Format_Grid.ColumnPrecision(1)=2
 ' ooi.Grid.Format_Grid.ColumnPrecision(3)=2
 ' ooi.Grid.Data_Grid.Text(1,1)="Absorb 310.8"
 ' ooi.Grid.Data_Grid.Text(1,2)="Absorb 313.1"
 ' ooi.Grid.Data_Grid.Text(1,3)="Absorb 315.4"

' For i = 1 To NWL1

' ooi.Grid.Data_Grid.value(1,i+3)=CWL1(i)

' Next i

' For i = 1 To NWL1 'NWL2

' ooi.Grid.Data_Grid.value(1,NWL1+3+i)=CWL2(i)

Resonance Ltd.	Operating Manual	Page 17 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

```
' Next i

' set numerical precision and format

'
' ooi.Grid.Format_Grid.ColumnPrecision(2)=2
' ooi.Grid.Format_Grid.ColumnPrecision(4)=2
' ooi.Grid.Format_Grid.ColumnExponentialFormat(2)=True
' ooi.Grid.Data_Grid.value(2,1)=1000000000
' ooi.Grid.Format_Grid.ColumnPrecision(4)=4
' ooi.Grid.Format_Grid.ColumnExponentialFormat(4)=True
' ooi.Grid.Data_Grid.value(2,1)=1000000000
' size the columns

' ooi.Grid.Format_Grid.AutosizeColumn(1,2,3,bfSizeBestFit)

'END comment out for excel

' initialize communication with excel

Set oXL=GetObject("Excel.Application")
Set oWB=oXL.ActiveWorkbook
Set oSheet=oWB.ActiveSheet

' put the wavelengths into the grid in the excel ss

oSheet.Cells(7,1)="Start Time"
oSheet.Cells(7,3)=Str(Time)
StartTimer = Timer
oSheet.Cells(10,1)="sec after start"
oSheet.Cells(10,2)="Absorb 310.8"
oSheet.Cells(10,3)="Absorb 313.1"
oSheet.Cells(10,4)="Absorb 315.4"

'
' For i = 1 To 237
'
' oSheet.Cells(10,18+i)=wl(Spect_Stpix+i*3)
'
' Next i
```

Resonance Ltd.	Operating Manual	Page 18 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

```
'      For i = 1 To NWL1
'
'          oSheet.Cells(10,4+i)=CWL1(i)
'
'      Next i

'
'      For i = 1 To NWL2
'
'          oSheet.Cells(10,4+NWL1+i)=CWL2(i)
'
'      Next i
```

```
'=====
'=====
'      DO NOT MODIFY THE FOLLOWING BLOCK OF CODE
'          Ready=INITDONE
'          END NONMODIFYABLE BLOCK
'=====
'=====
```

End Sub

Public Sub EventGenerator_NewSpectraReady(ByVal windowname As String, ByVal channels As Integer)

```
'=====
'=====
'      DO NOT MODIFY THE FOLLOWING BLOCK OF CODE
'          If Ready<>INITDONE Then Exit Sub
'          If ooi.InNextScanWait=True Then
'              ooi.InNextScanWait=False
'              Exit Sub
'          End If
'          If ooi.InTimeWait=True Then Exit Sub
'          If InEvent=True Then Exit Sub
'          InEvent=True
'          END NONMODIFYABLE BLOCK
'=====
'=====
```

```
'#####
'#####
```

Resonance Ltd.	Operating Manual	Page 19 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

```
' Insert your event code here
#####
#####
```

```
' this On Error statement is necessary. when you edit a cell
' in Excel while the script is running, Excel generates an
' error. automation clients are unalbe to access cell contents
' while it is being edited. the On Error statement causes the
' script to gracefully continue
```

On Error Resume Next

```
' with each new spectrum, get the desired wavelength
' Dim wl As Double
' wl=340 ' oSheet.Cells(5,3)
```

' define local variables

```
Dim M As Double
Dim X As Double
Dim B As Double
Dim RFlux1(9)
Dim CRFlux1(9)
Dim CRFlux2(5)
Dim RFlux2(5)
Dim CurPix As Integer
```

' get processed spectral data

```
ooi.SpectralData.GetProcessedSpectrum(chMaster,data(0))
```

***** Create compressed spectrum *****

'Create average

' integrate

For i =1 To NWL1

Resonance Ltd.	Operating Manual	Page 20 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

```

RFlux1(i) = 0
For CurPix = Stpix1(i) To Enpix1(i)
    RFlux1(i) = RFlux1(i) + data(CurPix)
Next CurPix
RFlux1(i) = RFlux1(i)/(Enpix1(i)-Stpix1(i))
Next i

```

```

For i = 1 To NWL2
    RFlux2(i) = 0
    For CurPix = Stpix2(i) To Enpix2(i)
        RFlux2(i) = RFlux2(i) + data(CurPix)
    Next CurPix
    RFlux2(i) = RFlux2(i)/(Enpix2(i)-Stpix2(i))
Next i

```

' USE CODE BELOW IF DATA GOING TO EXCEL

' put the data into the excel spreadsheet

j = j+1

M = (RFlux1(2)-RFlux1(1))/(CWL1(2)-CWL1(1))
B = RFlux1(1)

For i = 1 To 2

CRFlux1(i) = RFlux1(i)

Next i

For i = 3 To NWL1

X = CWL1(i) - CWL1(2)

CRFlux1(i) = RFlux1(i) - B *(M*X+ B)

Next i

Resonance Ltd.	Operating Manual	Page 21 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

For i = 1 To NWL2

$$X = \text{CWL2}(i) - \text{CWL1}(2)$$

$$\text{CRFlux2}(i) = \text{RFlux2}(i) - B'(M*X + B)$$

Next i

ooi.oSheet.Cells(10+j,1)= Str(Timer-StartTimer)

ooi.oSheet.Cells(10+j,2)= -

Log(2*CRFlux1(4)/(CRFlux1(3)+CRFlux1(5))/BaseRat309)

ooi.oSheet.Cells(10+j,3)= -

Log(2*CRFlux1(6)/(CRFlux1(5)+CRFlux1(7))/BaseRat311)

ooi.oSheet.Cells(10+j,4)= -

Log(2*CRFlux1(8)/(CRFlux1(7)+CRFlux1(9))/BaseRat313)

' For i = 1 To NWL1

' ooi.oSheet.Cells(10+j,4+i)= RFlux1(i)

' Next i

' For i = 1 To NWL2

' ooi.oSheet.Cells(10+j,i+4+NWL1)= RFlux2(i)

' Next i

'Put continuous spectrum into Excel

' For i = 1 To 237

' oSheet.Cells(10+j,18+i)=(data(Spect_Stpix+i*3-1)+data(Spect_Stpix+i*3)+data(Spect_Stpix+i*3+1))/3

' Next i

' USE CODE BELOW IF DATA GOING TO SHEET IN OOI32 PLAT

' put the data into the grid in the ooibase sheet

Resonance Ltd.	Operating Manual	Page 22 OF 22
Electro-Optics Research, Development and Manufacturing	Resonance Micro-DOAS RMDI	RMDI system.doc

```

'      ooi.Grid.Data_Grid.value(2,1) = -Log(2*RFlux1(4)/(RFlux1(3) +
RFlux1(5))/BaseRat309)
'      ooi.Grid.Data_Grid.value(2,2) = -Log(2*RFlux1(6)/(RFlux1(5) +
RFlux1(7))/BaseRat311)
'      ooi.Grid.Data_Grid.value(2,3) = -Log(2*RFlux1(8)/(RFlux1(7) +
RFlux1(9))/BaseRat313)
'      ooi.Grid.Data_Grid.value(3,1) = -
Log(2*CRFlux1(4)/(CRFlux1(3)+CRFlux1(5))/BaseRat309)
'      ooi.Grid.Data_Grid.value(3,2) = -
Log(2*CRFlux1(6)/(CRFlux1(5)+CRFlux1(7))/BaseRat311)
'      ooi.Grid.Data_Grid.value(3,3) = -
Log(2*CRFlux1(8)/(CRFlux1(7)+CRFlux1(9))/BaseRat313)

'      For i =1 To NWL1

'          ooi.Grid.Data_Grid.value(2,i+3)= RFlux1(i)
'          ooi.Grid.Data_Grid.value(3,i+3)= CRFlux1(i)

'      Next i

'      For i =1 To NWL2

'          ooi.Grid.Data_Grid.value(2,i+NWL1+3)= RFlux2(i)
'          ooi.Grid.Data_Grid.value(3,i+NWL1+3)= CRFlux2(i)

'      Next i

'=====
'=====
'      DO NOT MODIFY THE FOLLOWING BLOCK OF CODE
'          InEvent=False
'          END NONMODIFYABLE BLOCK
'=====
'=====
End Sub

```